



Chicken egg albumin (DAG2544)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95% , as determined by SDS-PAGE
---------------	-----------------------------------

Conjugate	Unconjugated
------------------	--------------

Applications	WB/ELISA
---------------------	----------

Molecular Weight	45 kDa
-------------------------	--------

Format	Lyophilized Powder
---------------	--------------------

Size	1 g
-------------	-----

Storage	Stored at 4°C
----------------	---------------

BACKGROUND

Introduction	Ovalbumin (abbreviated OVA) is the main protein found in egg white, making up 60-65% of the total protein. Ovalbumin displays sequence and three-dimensional homology to the serpin superfamily, but unlike most serpins it is not a serine protease inhibitor. The function of ovalbumin is unknown, although it is presumed to be a storage protein.
---------------------	--

Keywords	OVA; Ovalbumin
-----------------	----------------