



Chicken Trypsin Inhibitor (DAG2546)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Recombinant Trypsin Inhibitor Protein
Species	Chicken
Purity	Protein
Conjugate	Unconjugated
Format	Dried powder.
Concentration	Please see the vial label for concentration
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	The soybean trypsin inhibitor was first crystallized by Kunitz in 1945 and is one of several such inhibitors found in soybeans. It inhibits the proteolytic action approximately of an equal weight of crystalline trypsin by combining with trypsin to form a stable compound. Chymotrypsin is only slightly inhibited by soy inhibitor.
Keywords	25 kDa trypsin inhibitor (Protease inhibitor 15 preproprotein); CRISP8; KTi; Kti S; KTi3; KtiS; Kunitz trypsin inhibitor; Kunitz type trypsin inhibitor A; P24TI; P25TI; Trypsin inhibitor A