



Recombinant Human Human La/SSB Antigen [His] (DAG264)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	≥90% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA, CLIA, FCM
Molecular Weight	39 kDa
Reconstitution	Centrifuge before opening. Reconstitute with sterile distilled water to a concentration of 0.1-1.0 mg/mL. Aliquot and store at -80°C.
Format	Lyophilized
Size	100 µg, 1 mg
Buffer	Lyophilized from a 0.22 µm filtered PBS solution with 5% trehalose added as a protectant prior to lyophilization
Preservative	None
Storage	Store at -80°C

BACKGROUND

Keywords La/SSB; SSB (La); Sjögren syndrome type B antigen