



# Recombinant Human PAPP-A Protein [His] (DAG268)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Purified recombinant human Pregnancy Associated Plasma Protein-1 protein with an N-terminal hexahistidine tag.
<b>Species</b>	Human
<b>Conjugate</b>	His
<b>Applications</b>	Control, WB
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	20mM tris HCl with 50% glycerol
<b>Preservative</b>	None
<b>Storage</b>	Store at -20°C. Avoid Freeze/Thaw Cycles

## BACKGROUND

<b>Introduction</b>	Pregnancy-associated plasma protein A, pappalysin 1, also known as PAPP-A, is a protein used in screening tests for Down syndrome. This gene encodes a secreted metalloproteinase which cleaves insulin-like growth factor binding proteins (IGFBPs).
<b>Keywords</b>	PAPP-A1; PAPP-A; Pregnancy Associated Plasma Protein-1; Pregnancy Associated Plasma Protein; Pregnancy-associated plasma protein A