



EDDP [BTG] (DAG2983)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	EDDP BTG conjugate
Species	N/A
Conjugate	BTG
Applications	immunoassay development or other applications.
Format	Liquid
Concentration	4.0 mg/mL
Size	1 mg
Buffer	Supplied in 20 mM potassium phosphate, pH 7.5
Preservative	0.1% Sodium Azide
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	It is the main metabolite of methadone.
Keywords	BAY 78418; bay78418; BAYER 78418; bayer78418; Blastoff; Dithiophosphorsaeure-O-aethyl-S,S-diphenylester; Ediphenophos; Ediphenphos; Lutrol; O-Aethyl-S,S-diphenyl-dithiophosphat; Phosphorodithioic acid, O-ethyl S,S-diphenyl ester; phosphorodithioicacid