



## Methylphenidate [BSA] (DAG3007)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Product Overview	Methylphenidate, BSA-Conjugated
Species	N/A
Conjugate	BSA
Applications	immunoassay development or other applications.
Format	Liquid
Size	1 mg, 5 mg
Buffer	Conjugate supplied in PBS, pH 7.5
Preservative	0.01% Sodium Azide
Storage	2-8°C short term, -20°C long term

### BACKGROUND

Introduction	Methylphenidate is a central nervous system stimulant. It affects chemicals in the brain and nerves that contribute to hyperactivity and impulse control. Methylphenidate is used to treat attention deficit disorder (ADD) and attention deficit hyperactivity disorder (ADHD). It is also used to treat uncontrollable periods of daytime sleep (narcolepsy).
Keywords	MPD; MPH; Ritalin; Methylphenidate