



Nortriptyline [BSA] (DAG3012)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Nortriptyline, BSA-Conjugated
Species	N/A
Conjugate	BSA
Applications	immunoassay development or other applications.
Format	Liquid
Concentration	1.0 mg/mL
Size	1 mg
Buffer	Supplied in 0.015 M phosphate, 0.15 M NaCl, pH 7.2
Preservative	0.1% Sodium Azide
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	Nortriptyline is a second-generation tricyclic antidepressant (TCA).
Keywords	Nortriptyline; Sensoval; Aventyl; Pamelor; Norpress; Allegron; Noritren; Nortrilen