



## Progesterone-11 [HRP] (DAG3019H)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Progesterone also known as P4 (pregn-4-ene-3,20-dione) is a C-21 steroid hormone involved in the female menstrual cycle, pregnancy (supports gestation) and embryogenesis of humans and other species. Progesterone belongs to a class of hormones called progestogens, and is the major naturally occurring human progestogen.
<b>Species</b>	N/A
<b>Purity</b>	Highly purified by gel filtration
<b>Conjugate</b>	HRP
<b>Applications</b>	ELISA ( $\geq 1:1000$ )
<b>Format</b>	Liquid
<b>Concentration</b>	Please refer to the vial label for the specific concentration.
<b>Size</b>	1 ml
<b>Buffer</b>	0.1M PBS with 0.1% Tween-20, 0.5% solubilized casein and dye
<b>Storage</b>	2-8°C. Do not freeze.
<b>Ship</b>	Wet ice