



Testosterone-19 [BSA] (DAG3024)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Testosterone is a steroid hormone from the androgen group and is found in mammals, reptiles, birds, and other vertebrates. In mammals, testosterone is secreted primarily by the testicles of males and the ovaries of females, although small amounts are also secreted by the adrenal glands. It is the principal male sex hormone and an anabolic steroid.
Species	N/A
Conjugate	BSA
Applications	N/A
Format	Liquid
Concentration	≥ 1.0 mg/mL
Size	1 mg
Buffer	Phosphate buffered saline, pH 7.2
Preservative	0.1% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.

BACKGROUND

Keywords Testosterone, Testosterone-19