



Theophylline [HRP] (DAG3030)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Theophylline, HRP-Conjugated
Species	N/A
Conjugate	HRP
Applications	immunoassay development or other applications.
Format	Liquid
Concentration	1 mg/mL
Size	1 ml
Buffer	1:1 glycerol / 0.015 M phosphate, 0.15 M NaCl, pH 7.2
Preservative	None
Storage	Store at -20°C.

BACKGROUND

Introduction	Theophylline is a methylxanthine drug used in therapy for respiratory diseases such as COPD and asthma under a variety of brand names.
Keywords	Theophylline; 1 3-dimethylxanthine; Theolair; Slo-Bid