



Triiodothyronine (T3) [BSA] (DAG3033)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Triiodothyronine, BSA-Conjugated
Species	N/A
Conjugate	BSA
Applications	immunoassay development or other applications.
Reconstitution	Reconstitute with laboratory grade water.
Format	Lyophilized
Size	1mg, 25mg
Preservative	None
Storage	Aliquot and store at -20 °C. Avoid repeated Freeze/Thaw cycles

BACKGROUND

Introduction	The thyroid hormones, thyroxine (T4) and triiodothyronine (T3), are tyrosine-based hormones produced by the thyroid gland primarily responsible for regulation of metabolism. An important component in the synthesis of thyroid hormones is iodine. The major form of thyroid hormone in the blood is thyroxine (T4), which has a longer half life than T3.
--------------	--