



Digoxin [BSA] (DAG3043)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Digitalis lanata Digoxin, BSA-Conjugated
Species	N/A
Conjugate	BSA
Applications	immunoassay development or other applications.
Format	Liquid
Concentration	1.0 mg/mL
Size	1 mg
Buffer	Supplied in 0.015 M phosphate, 0.15 M NaCl, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Store at 2-8°C short term and -20°C long term

BACKGROUND

Introduction	Digoxin is a beta blocker drug originally derived from the foxglove plant, Digitalis lanata. Digoxin is used primarily to improve the pumping ability of the heart in congestive heart failure (CHF), and treat problems such as high blood pressure. It is also used to help normalize some dysrhythmias (abnormal types of heartbeat).
Keywords	Digitalis; Lanoxin; Digoxin