



Recombinant Protein A [HIS] (DAG3055)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Recombinant VanA fusion protein, HIS-tagged
Species	N/A
Purity	> 95% by SDS-PAGE
Conjugate	HIS
Applications	immunoassay development or other applications.
Buffer	TBS, pH 8.0 with 0.09% Sodium Azide added as preservative
Preservative	0.09% Sodium Azide
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	VanA (Vancomycin Resistant Protein A) of Vancomycin resistant Enterococcus sp. (VRE)
Keywords	Vancomycin Resistant Protein A; VanA