



Native CMV Infected Cell Extract (Strain AD169) (DAG3083)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	N/A
Applications	ELISA, WB
Format	Liquid
Size	1 mL, 25 mL
Buffer	100mM Glycine, pH 9.5
Preservative	None
Storage	Store at -80°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Keywords	Cytomegalovirus; CMV; HCMV; HHV-5
----------	-----------------------------------