



Secobarbital [BSA] (DAG3102)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Secobarbital, BSA-Conjugated
Species	N/A
Conjugate	BSA
Applications	Lateral Flow
Format	Liquid
Concentration	1.21mg/mL
Size	1 mg
Buffer	Supplied in distilled water
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	Secobarbital is a barbiturate derivative. Barbiturates are a class of around 12 compounds derivatised from barbituric acid. They are central nervous system (CNS) depressants and can be used as sedatives, hypnotics, anaesthetics and anti-epileptic drugs.
Keywords	Secobarbital