



# Tobramycin [BSA] (DAG3141)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	Tobramycin, BSA-Conjugated
Specificity	14/1 (Moles Antigen/Moles Protein)
Species	N/A
Conjugate	BSA
Applications	immunogen
Preservative	None
Storage	2-8°C short term, -20°C long term

## BACKGROUND

Introduction	Tobramycin is an aminoglycoside antibiotic used to treat infections caused by many different bacteria. Tobramycin is usually administered by intravenous infusion, intramuscular injection, or inhalation.
Keywords	Tobramycin