



Native CMV (Strain AD169) Antigen (DAG3225)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	CMV
Purity	Density gradient purified
Conjugate	Unconjugated
Applications	IgM Assays: ELISA, WB
Format	Liquid
Size	100 µg
Buffer	10 mM Tris, 100 mM NaCl, 1 mM EDTA, 0.1% Triton X-100, pH 8
Preservative	None
Storage	Store at -70°C or below.

BACKGROUND

Introduction	Cytomegalovirus is a member of the herpes virus group, which includes herpes simplex virus types 1 and 2, varicella zoster virus (which causes chicken pox), and Epstein–Barr virus (which causes infectious mononucleosis). These viruses share a characteristic ability to remain dormant within the body over a long period.
Keywords	Cytomegalovirus; CMV