



# 1,3-Butadiene [BSA] (DAG3250)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	1,3-Butadiene, BSA-conjugated
Specificity	1,3-Butadiene conjugated with bovine serum albumin (BSA).
Species	N/A
Conjugate	BSA
Applications	immunohistochemistry and immunocytochemistry.
Reconstitution	Reconstitute with deionized water (250 µg)
Format	Lyophilized
Concentration	Not applicable
Size	1 mg
Preservative	None
Storage	2-8°C short term, -20°C long term

## BACKGROUND

Introduction	1,3-Butadiene is a simple conjugated diene. It is an important industrial chemical used as a monomer in the production of synthetic rubber.
Keywords	1,3-Butadiene