



3-Hydroxy-Kyurenine [BSA] (DAG3251)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	3-Hydroxy-Kyurenine (BSA) conjugate
Specificity	3-Hydroxykynurenine conjugated to bovine serum albumin (BSA).
Species	N/A
Conjugate	BSA
Applications	immunohistochemistry and immunocytochemistry.
Reconstitution	Reconstitute with deionized water (250 µg)
Format	Lyophilized
Concentration	Not applicable
Size	1 mg
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Keywords	3-Hydroxy-Kyurenine
----------	---------------------