



3-Indolacetic Acid [BSA] (DAG3252)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	3-Indolacetic Acid, BSA-conjugated
Specificity	3-Indolacetic Acid conjugated with bovine serum albumin (BSA).
Species	N/A
Conjugate	BSA
Applications	immunohistochemistry and immunocytochemistry.
Reconstitution	Reconstitute with deionized water (250 µg)
Format	Lyophilized
Concentration	Not applicable
Size	1 mg
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Keywords	3-Indolacetic Acid
----------	--------------------