



4,5-Dihydroxytryptamine [BSA] (DAG3253)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	4,5-Dihydroxytryptamine, BSA-conjugated
Specificity	4,5-Dihydroxytryptamine conjugated with glutaraldehyde (G) and bovine serum albumin (BSA).
Species	N/A
Conjugate	BSA
Applications	immunohistochemistry and immunocytochemistry.
Reconstitution	Reconstitute with deionized water (250 µg)
Format	Lyophilized
Concentration	Not applicable
Size	1 mg
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	Tryptamine substituted with two hydroxyl groups in positions 4 and 5.
Keywords	4,5-DHT; 4,5-dihydroxytryptamine