



5-Hydroxy-Indole Acetic Acid [BSA] (DAG3258)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	5-Hydroxy-Indole Acetic Acid, BSA-conjugated
Specificity	5-Hydroxy-Indole Acetic Acid (5HIAA) conjugated with bovine serum albumin (BSA).
Species	N/A
Conjugate	BSA
Applications	immunohistochemistry and immunocytochemistry
Reconstitution	Reconstitute with deionized water (250 µg)
Format	Lyophilized
Concentration	Not applicable
Size	1 mg
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Keywords	5-Hydroxy-Indole
----------	------------------