



Kainic Acid [BSA] (DAG3357)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10 mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C

BACKGROUND

Keywords	Kainic acid; Kainate; (3S,4S)-3-(Carboxymethyl)-4-(prop-1-en-2-yl)-L-proline
----------	--