



L-Kynurenine [BSA] (DAG3359)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10 mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C

BACKGROUND

Keywords	L-Kynurenine; (S)-Kynurenine; (2S)-2-Amino-4-(2-aminophenyl)-4-oxo-butanoic acid
----------	--