



## Melatonin [AG-BSA] (DAG3367)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|              |                   |
|--------------|-------------------|
| Conjugate    | BSA               |
| Applications | ELISA, LFIA       |
| Format       | Liquid            |
| Size         | 1 mg              |
| Buffer       | 10 mM PBS, pH 7.4 |
| Preservative | None              |
| Storage      | Store at -20°C    |

### BACKGROUND

|          |                                                          |
|----------|----------------------------------------------------------|
| Keywords | Melatonin; N-[2-(5-methoxy-1H-indol-3-yl)ethyl]acetamide |
|----------|----------------------------------------------------------|