



Myristic Acid [BSA] (DAG3369)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Myristic Acid, BSA-conjugated
Specificity	Myristic Acid conjugated with bovine serum albumin (BSA).
Species	N/A
Conjugate	BSA
Applications	immunohistochemistry and immunocytochemistry
Reconstitution	Reconstituted in deionized water (250 µg)
Format	Lyophilized
Size	1 mg
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	Myristic acid, also called tetradecanoic acid, is a common saturated fatty acid with the molecular formula $\text{CH}_3(\text{CH}_2)_{12}\text{COOH}$. A myristate is a salt or ester of myristic acid. Myristic acid is named after the nutmeg <i>Myristica fragrans</i> .
Keywords	Myristic acid