



NO₂-Tyrosine [G-BSA] (DAG3372)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	NO ₂ -Tyrosine, G-BSA-conjugated
Specificity	NO ₂ -Tyrosine conjugated with glutaraldehyde (G) and bovine serum albumin (BSA).
Species	N/A
Conjugate	G-BSA
Applications	immunohistochemistry and immunocytochemistry
Reconstitution	Reconstituted in deionized water (250 µg)
Format	Lyophilized
Size	1 mg
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Keywords	NO ₂ -Tyrosine
-----------------	---------------------------