



## NO-Dopamine [G-BSA] (DAG3376)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Product Overview	NO-Dopamine, G-BSA-conjugated
Specificity	NO-Dopamine conjugated with glutaraldehyde (G) and bovine serum albumin (BSA).
Species	N/A
Purity	Purity is > 95% pure as determined by 10% SDS-PAGE followed by coomassie staining
Conjugate	G-BSA
Applications	immunohistochemistry and immunocytochemistry
Reconstitution	Reconstituted in deionized water (250 µg)
Format	Lyophilized
Size	1 mg
Preservative	None
Storage	2-8°C short term, -20°C long term

### BACKGROUND

Keywords	Dopamine
----------	----------