



## NO-L-Glutamine [G-BSA] (DAG3381)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Product Overview	NO-L-Glutamine, G-BSA-conjugated
Specificity	NO-L-Glutamine conjugated with glutaraldehyde (G) and bovine serum albumin (BSA), NO = nitrosylated.
Species	N/A
Conjugate	G-BSA
Applications	immunohistochemistry and immunocytochemistry
Reconstitution	Reconstituted in deionized water (250 µg)
Format	Lyophilized
Size	1 mg
Preservative	None
Storage	2-8°C short term, -20°C long term

### BACKGROUND

Keywords	NO-L-Glutamine
----------	----------------