



Valproic Acid [ALP] (DAG3422A)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Valproic acid (VPA, valproate), an acidic chemical compound, has found clinical use as an anticonvulsant and mood-stabilizing drug, primarily in the treatment of epilepsy, bipolar disorder and prevention of migraine headaches.
Species	N/A
Conjugate	ALP
Applications	Suitable for use in immunoassay development or other applications.
Format	Liquid
Concentration	0.5 mg/mL
Size	1 mg
Buffer	0.02 M Tris HCl, pH 8.0
Storage	For long-term storage, aliquot and store at -20°C
Ship	Wet ice

BACKGROUND

Keywords	VPA; Valproic acid
----------	--------------------