



Chicken Insulin Like Growth Factor II (DAG3433)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95% as determined by HPLC
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Lyophilized white powder
Concentration	Batch dependent - please inquire should you have specific requirements
Size	100 µg
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	This gene is a member of the leukocyte immunoglobulin-like receptor (LIR) family, which is found in a gene cluster at chromosomal region 19q13.4. The encoded protein belongs to the subfamily B class of LIR receptors which contain two or four extracellular
Keywords	IGF2R; Insulin-like growth factor 2 receptor; MPR1