



Methyltetrahydrophthalic anhydride [BSA] (DAG343S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Methyltetrahydrophthalic anhydride, covalently linked to bovine serum albumin
Species	N/A
Purity	chromatography, dialysis
Conjugate	BSA
Applications	e.g. allergy test, immunoassays etc.
Molecular Weight	app. 65 kDa
Format	Liquid
Concentration	1 mg/mL, conjugation ratio is app. 1:10 - 1:40.
Size	1 mg
Buffer	PBS
Preservative	0.05% sodium azide
Storage	Store at 4°C.

BACKGROUND

Introduction	A cyclic dicarboxylic anhydride that is the cyclic anhydride of methyltetrahydrophthalic acid.
Keywords	drug-BSAconjugate, drug-conjugate, Methyltetrahydrophthalic anhydride BSAconjugate