



Human Vascular Endothelial Growth Factor A Polypeptid (DAG348)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Conjugate	Unconjugated
Molecular Weight	38.2kDa
Format	Lyophilized
Size	1 mg
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	Vascular endothelial growth factor A (VEGF-A) is a protein that in humans is encoded by the VEGFA gene.
Keywords	VEGFA; vascular endothelial growth factor A; VPF; VEGF; MVCD1; VEGF-A; MGC70609