



Recombinant HCV Nonstructural Protein 5 (a.a. 2061-2302) [Biotin] (DAG3503)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Recombinant Hepatitis C Virus NS5 (a.a. 2061-2302), Biotin-Conjugated
Specificity	Immunoreactive with sera of HCV-infected individuals.
Species	HCV
Purity	HCV-NS5 protein is > 95% pure as determined by 10% PAGE (coomassie staining).
Conjugate	Biotin
Applications	HCV-NS5 antigen is suitable for ELISA and Western blots, excellent antigen for detection of HCV with minimal specificity problems.
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	Hepatitis C Virus is a positive, single stranded RNA virus in the Flaviviridae family. The genome is approximately 10, 000 nucleotides and encodes a single polypeptide of about 3, 000 amino acids. The polypeptide is processed by host cell and viral proteases into three major structural proteins and several non structural proteins necessary for viral replication. Several different genotypes of HCV with slightly different genomic sequences have since been identified that correlate with differences in response to treatment with interferon alpha.
Keywords	HCV NS5; Hepatitis C virus non structural protein 5; Hepatitis C virus nonstructural protein 5; Hepatitis C Virus NS5; HCV NS5 Genotype 2a; Hepatitis C Virus NS5 Genotype 2a