



Recombinant *P. auratus* Insulin-like growth factor 1 receptor (DAG3577)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Recombinant Bream Insulin-like growth factor 1 receptor Grade
Species	<i>P. auratus</i>
Purity	> 90% (by HPLC and N-terminal sequence analysis)
Conjugate	Unconjugated
Format	Lyophilized white powder Dried from 0.1 M acetic acid and stored under nitrogen at a slight vacuum
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	This receptor binds insulin-like growth factor with a high affinity. It has tyrosine kinase activity. The insulin-like growth factor I receptor plays a critical role in transformation events. Cleavage of the precursor generates alpha and beta subunits. It is highly overexpressed in most malignant tissues where it functions as an anti-apoptotic agent by enhancing cell survival. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
Keywords	IGF1R; insulin-like growth factor 1 receptor; CD221; IGFR; IGFR; JTK13; MGC18216; IGF-I receptor; soluble IGF1R variant 1; soluble IGF1R variant 2; insulin-like growth factor I receptor; MGC142170; MGC142172; dalotuzumab; 1005389-60-5