



## Enoxazin [BSA] (DAG362S)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|                     |                                                     |
|---------------------|-----------------------------------------------------|
| Antigen Description | Enoxazin, covalently linked to bovine serum albumin |
| Species             | N/A                                                 |
| Purity              | chromatography, dialysis                            |
| Conjugate           | BSA                                                 |
| Applications        | e.g. allergy test, immunoassays etc.                |
| Molecular Weight    | app. 65 kDa                                         |
| Format              | Liquid                                              |
| Concentration       | 1 mg/mL, conjugation ratio is app. 1:10 - 1:40.     |
| Size                | 1 mg                                                |
| Buffer              | PBS                                                 |
| Preservative        | 0.05% sodium azide                                  |
| Storage             | Store at 4°C.                                       |

### BACKGROUND

|          |                                                            |
|----------|------------------------------------------------------------|
| Keywords | drug-BSA conjugate, drug-conjugate, Enoxazin BSA conjugate |
|----------|------------------------------------------------------------|