



Caprylic acid [BSA] (DAG3646)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|------------------|---|
| Product Overview | Caprylic acid, BSA-Conjugated |
| Species | N/A |
| Conjugate | BSA |
| Applications | Immunohistochemistry, immunocytochemistry, ELISA, The optimal concentration should be determined by the user for each specific application. |
| Preservative | None |
| Storage | 2-8°C short term, -20°C long term |

BACKGROUND

| | |
|--------------|---|
| Introduction | Caprylic acid is the common name for the eight-carbon saturated fatty acid known by the systematic name octanoic acid. Its compounds are found naturally in the milk of various mammals, and as a minor constituent of coconut oil and palm kernel oil. It is an oily liquid that is minimally soluble in water with a slightly unpleasant rancid-like smell and taste. |
| Keywords | Caprylic acid |