



Azithromycin [BSA] (DAG382S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Azithromycin, covalently linked to bovine serum albumin
Species	N/A
Purity	chromatography, dialysis
Conjugate	BSA
Applications	e.g. allergy test, immunoassays etc.
Molecular Weight	app. 65 kDa
Format	Liquid
Concentration	1 mg/mL, conjugation ratio is app. 1:10 - 1:40.
Size	1 mg
Buffer	PBS
Preservative	0.05% sodium azide
Storage	Store at 4°C.

BACKGROUND

Introduction	Azithromycin is an azalide, derived from erythromycin, and a member of a subclass of macrolide antibiotics with bacteriocidal and bacteriostatic activities. Azithromycin reversibly binds to the 50S ribosomal subunit of the 70S ribosome of sensitive microorganisms, thereby inhibiting the translocation step of protein synthesis, wherein a newly synthesized peptidyl tRNA molecule moves from the acceptor site on the ribosome to the peptidyl (donor) site, and consequently inhibiting RNA-dependent protein synthesis leading to cell growth inhibition and
---------------------	--

cell death.

Keywords

drug-BSA conjugate, drug-conjugate, Azithromycin BSA conjugate
