



Propoxyphene [BSA] (DAG397)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------|--|
| Conjugate | BSA |
| Applications | N/A |
| Format | Liquid |
| Concentration | Batch dependent - please inquire should you have specific requirements |
| Size | 1 mg |
| Buffer | Distilled water |
| Preservative | None |
| Storage | 2-8°C short term, -20°C long term |

BACKGROUND

| | |
|--------------|---|
| Introduction | Dextropropoxyphene is an analgesic in the opioid category, patented in 1955 and manufactured by Eli Lilly and Company. It is an optical isomer of levopropoxyphene. It is intended to treat mild pain and also has antitussive (cough suppressant) and local anaesthetic effects. The drug has been taken off the market in Europe and the US due to concerns of fatal overdoses and heart arrhythmias. It is still available in Australia, albeit with restrictions after an application by its manufacturer to review its proposed banning. |
| Keywords | Dextropropoxyphene; Analgesia |