



Recombinant HCV type 6a Nonstructural Protein 5 [GST] (DAG4007)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Recombinant Hepatitis C Virus Non-structural Protein 5 Genotype-6a, GST-Tagged
Species	HCV
Purity	greater than 95% by 10% SDS-PAGE
Conjugate	GST
Applications	ELISA and Western blots, excellent antigen for detection of HCV with minimal specificity problems.
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	Hepatitis C Virus is a positive, single stranded RNA virus in the Flaviviridae family. The genome is approximately 10, 000 nucleotides and encodes a single polyprotein of about 3, 000 amino acids. The polyprotein is processed by host cell and viral proteases into three major structural proteins and several non structural proteins necessary for viral replication. Several different genotypes of HCV with slightly different genomic sequences have since been identified that correlate with differences in response to treatment with interferon alpha.
Keywords	HCV NS5; Hepatitis C virus non structural protein 5; Hepatitis C virus nonstructural protein 5; Hepatitis C Virus NS5; HCV NS5 Genotype 2a; Hepatitis C Virus NS5 Genotype 2a