



# Recombinant HEV Mosaic Protein (HEV-1+HEV-4) (DAG4011)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Conjugate	TBD
Applications	ELISA, CLIA, LFIA
Molecular Weight	41.9KD, 24.5KD
Reconstitution	Reconstitute at 1mg/mL in sterile water
Format	Lyophilized
Size	1 mg
Buffer	25mM Tris-HCl, 0.02%SDS, pH8.5
Preservative	None
Storage	Store at 2°C to 8°C short term or frozen below -20°C

## BACKGROUND

Keywords	HEV; Hepatitis E virus; HEV-1; HEV-4
----------	--------------------------------------