



## Methadone [BSA] (DAG410)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|                  |   |
|------------------|---|
| Product Overview | Methadone-BovineSerum Albumin (MTD-BSA) Conjugate               |
| Species          | N/A   |
| Conjugate        | BSA   |
| Applications     | Specific methodologies have not been tested using this product. |
| Format           | Liquid  |
| Concentration    | 5.1 mg/mL   |
| Size             | 1 mg  |
| Buffer           | PBS, pH 7.2   |
| Preservative     | 0.02% Sodium Azide  |
| Storage          | 2-8°C short term, -20°C long term                               |

### BACKGROUND

|              |   |
|--------------|---|
| Introduction | Methadone is a synthetic opioid, used medically as an analgesic and a maintenance anti-addictive for use in patients on opioids. Although chemically unlike morphine or heroin, methadone also acts on the opioid receptors and thus produces many of the same effects. |
| Keywords     | Methadone   |