



Penicillin V [HRP] (DAG414S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Penicillin V, covalently linked to horseradish peroxidase
Species	N/A
Purity	chromatography, dialysis
Conjugate	HRP
Applications	e.g. allergy test, immunoassays etc.
Molecular Weight	app. 44 kDa
Reconstitution	Reconstitute with 1 ml PBS without azide.
Format	Lyophilized
Concentration	1 mg/mL working concentration
Size	1 mg
Preservative	None
Storage	Store at -20°C.

BACKGROUND

Introduction	Penicillin is a group of Beta-lactam antibiotics used in the treatment of bacterial infections caused by susceptible, usually Gram-positive, organisms. beta-lactam antibiotics work by inhibiting the formation of peptidoglycan cross-links in the bacterial cell wall, which results in cytolysis.
Keywords	drug-HRP conjugate, drug-conjugate, penicillin V HRP conjugate