



Metoprolol [HRP] (DAG423S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Metoprolol, covalently linked to horseradish peroxidase
Species	N/A
Purity	chromatography, dialysis
Conjugate	HRP
Applications	e.g. allergy test, immunoassays etc.
Molecular Weight	app. 44 kDa
Reconstitution	Reconstitute with 1 ml PBS without azide.
Format	Lyophilized
Concentration	1 mg/mL working concentration
Size	1 mg
Preservative	None
Storage	Store at -20°C.

BACKGROUND

Introduction	Metoprolol is a cardioselective competitive beta-1 adrenergic receptor antagonist with antihypertensive properties and devoid of intrinsic sympathomimetic activity. Metoprolol antagonizes beta 1-adrenergic receptors in the myocardium, thereby reducing the rate and force of myocardial contraction leading to a reduction in cardiac output. This agent may also reduce the secretion of renin with subsequent reduction in levels of angiotensin II thereby preventing vasoconstriction and aldosterone secretion.
--------------	---

Keywords

drug-HRP conjugate, drug-conjugate, Metoprolol HRP conjugate
