



Theophylline [BSA] (DAG426)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Theophylline-Bovine Serum Albumin Conjugate
Species	N/A
Conjugate	BSA
Applications	Specific methodologies have not been tested using this product.
Reconstitution	Reconstitute with 1.5ml laboratory grade water.
Format	Lyophilized.
Concentration	Not applicable
Size	1 mg
Buffer	Lyophilized from 7mM PBS
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	Theophylline is a methylxanthine drug used in therapy for respiratory diseases such as COPD and asthma under a variety of brand names.
Keywords	Theophylline; 1,3-dimethylxanthine; Theolair; Slo-Bid