



Diatrizoate [HRP] (DAG436S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|----------------------------|--|
| Antigen Description | Diatrizoate, covalently linked to horseradish peroxidase |
|----------------------------|--|

| | |
|----------------|-----|
| Species | N/A |
|----------------|-----|

| | |
|---------------|--------------------------|
| Purity | chromatography, dialysis |
|---------------|--------------------------|

| | |
|------------------|-----|
| Conjugate | HRP |
|------------------|-----|

| | |
|---------------------|--------------------------------------|
| Applications | e.g. allergy test, immunoassays etc. |
|---------------------|--------------------------------------|

| | |
|-------------------------|-------------|
| Molecular Weight | app. 44 kDa |
|-------------------------|-------------|

| | |
|-----------------------|---|
| Reconstitution | Reconstitute with 1 ml PBS without azide. |
|-----------------------|---|

| | |
|---------------|-------------|
| Format | Lyophilized |
|---------------|-------------|

| | |
|----------------------|-------------------------------|
| Concentration | 1 mg/mL working concentration |
|----------------------|-------------------------------|

| | |
|-------------|------|
| Size | 1 mg |
|-------------|------|

| | |
|---------------------|------|
| Preservative | None |
|---------------------|------|

| | |
|----------------|-----------------|
| Storage | Store at -20°C. |
|----------------|-----------------|

BACKGROUND

| | |
|---------------------|---|
| Introduction | A commonly used x-ray contrast medium. As DIATRIZOATE MEGLUMINE and as Diatrizoate sodium, it is used for gastrointestinal studies, angiography, and urography. |
|---------------------|---|

| | |
|-----------------|---|
| Keywords | drug-HRP conjugate, drug-conjugate, Diatrizoate HRP conjugate |
|-----------------|---|