



Recombinant *D. melanogaster* Importin α [His] (DAG4423)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	<i>Drosophila melanogaster</i> Importin α , His-tagged
Species	<i>D. melanogaster</i>
Purity	> 95% by SDS-PAGE
Conjugate	His
Format	Liquid. Supplied as a 1 mg/ml solution in 50 mM Tris-HCl, pH 7.5, 2 mM DTE, and 5 % Glycerol
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction Importin alpha, or karyopherin alpha proteins are involved in the import of proteins into the cell nucleus. They are a sub-family of karyopherin proteins.

Keywords Importin alpha; Importin α ; *Drosophila melanogaster*; *D. melanogaster*; *Drosophila*