



Recombinant HCMV (AB 32/32) [His] (DAG4580)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Recombinant Human cytomegalovirus AB 32/32, His-tagged
Species	HCMV
Conjugate	His
Buffer	100 mM KH ₂ PO ₄ , 100 mM Na ₂ HPO ₄ ·2H ₂ O, 300 mM NaCl, 10% glycerol, 0,05 mM mercaptoethanol, pH = 8
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	Cytomegalovirus (CMV), is a genus of Herpes viruses; in humans the species is known as Human herpesvirus 5 (HHV-5). It belongs to the Betaherpesvirinae subfamily of Herpesviridae. The name means "very big cell virus".
Keywords	Group I (dsDNA); Herpesviridae; Betaherpesvirinae; Cytomegalovirus; CMV; CMV glycoprotein B; CMV glycoprotein GP55; Cytomegalovirus glycoprotein GP55; gB; Glycoprotein B; Glycoprotein GP55; HCMV; HCMV glycoprotein B; HCMV glycoprotein GP55; HHV 5; HHV 5 g