



Chicken Leptin Binding Domain (DAG4593)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Recombinant Chicken Leptin Binding Domain
Species	Chicken
Conjugate	Unconjugated
Format	Sterile filtered, lyophilized
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	Leptin Binding Domain Chicken Recombinant also called Leptin Receptor produced in E.coli is a single, non-glycosylated, polypeptide chain containing 208 amino acids and having a molecular mass of 24.5 kDa. Chicken Leptin Receptor consists of the cytokine binding domain of leptin receptor amino acids 420-626 of chicken leptin receptor. The Leptin Binding Domain is purified by proprietary chromatographic techniques
Keywords	Chicken Leptin Binding Domain protein; Chicken; Leptin Binding Domain; OB Protein protein; Obesity Protein protein; OBS protein; Obesity factor protein; Leptin Receptor protein