



## Alcuronium [HRP] (DAG460S)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Antigen Description	Alcuronium, covalently linked to horseradish peroxidase
Species	N/A
Purity	chromatography, dialysis
Conjugate	HRP
Applications	e.g. allergy test, immunoassays etc.
Molecular Weight	app. 44 kDa
Reconstitution	Reconstitute with 1 ml PBS without azide.
Format	Lyophilized
Concentration	1 mg/mL working concentration
Size	1 mg
Preservative	None
Storage	Store at -20°C.

### BACKGROUND

Introduction	Alcuronium is an indole alkaloid of the curare family. A neuroblocker, it is often used in chloride form as an anesthesia adjuvant. It has a role as a muscle relaxant.
Keywords	drug-HRP conjugate, drug-conjugate, Alcuronium HRP conjugate